

OASYS® OTS

OAsys® is a double upright brace that provides medial or lateral compartment unloading for the relief of OA knee pain, plus functional support for patients with ligament instability. Equipped with Accutrac® ROM hinges¹ and an Acculign™ slide bar², it is also suitable for post-op and long-term rehab [>4weeks] protocols requiring unloading and/or adjustable flexion/extension control. Quick-release buckles³ make OAsys easy to apply and remove.



INDICATIONS

- Medial or lateral compartment unloading for unicompartmental osteoarthritis or articular cartilage healing
- Post-op, rehab and functional support plus unloading for protocols requiring reduced medial or lateral joint loads during recovery
- Functional support with or without varus/valgus alignment
- ACL, MCL, LCL, PCL, rotary and combined instabilities
- Low and medium contact/impact activities

PART#	LEG	SIZE	MEASURING ¹	COLOR
OA3SR	Right	Small	Caliper: 3.375-3.875"	Silver, Champagne
OA3SL	Left	Small	Circ: 11-12.25"	Clear Carbon
OA3MR	Right	Medium	Caliper: 3.875-4.375"	Silver, Champagne
OA3ML	Left	Medium	Circ: 12.25-13.875"	Clear Carbon
OA3LR	Right	Large	Caliper: 4.375-5"	Silver, Champagne
OA3LL	Left	Large	Circ: 13.875-15.875"	Clear Carbon
OA3XLR	Right	XLarge	Caliper: 5-5.5"	Silver, Champagne
OA3XLL	Left	XLarge	Circ: 15.875-17.375"	Clear Carbon

¹ Caliper: M/L at joint width in a standard, weight-bearing position
Circumference: just below the distal border of the patella



Medial View
showing
Acculign™
slide bar